

TED

TED Anti-embolism stockings

- ~ Colour: white
- ~ Inspection hole at toe
- ~ Sizing refer TED sizing charts
- ~ Order by product code





TED

TED knee length – sizing chart #1

Calf size	Calf (cm) circ	Leg length, heel to knee (cm)			
		Regular	Long		
small	up to 30.5	up to 40.6	40.6+		
medium	30.5-38.1	up to 43.2	43.2+		
large	38.1-44.5	up to 45.7	45.7+		
x large	44.5-50.8	up to 45.7	45.7+		

TED thigh length – sizing chart #2

Calf size	Calf circ (cm)	Thigh circ (cm)	Leg length, heel to groin (cm)		
			Short	Regular	Long
small	up to 30.5	up to 63.5	up to 73.7	73.7-83.8	83.8+
medium	30.5-38.1	up to 63.5	up to 73.7	73.7-83.8	83.8+
large	38.1-44.5	up to 63.5	up to 73.7	73.7-83.8	83.8+

TED thigh length with belt – sizing chart #3

Calf size	Calf circ (cm)	Thigh circ (cm)	Leg length, heel to groin (cm)	
			Regular	Long
x small	up to 24.5	up to 63.5	up to 71.1	71.1+
small	24.5-30.5	up to 63.5	up to 72.4	72.4+
medium	30.5-38.1	up to 63.5	up to 72.4	72.4+
large	38.1-44.5	up to 63.5	up to 73.7	73.7+
x large	38.1-44.5	63.5- 81.3	up to 72.4	72.4+